

**ANDITORK First – Bottle top/cap torque tester**

↳ ANDITORK is designed for measuring the opening/closing torque of bottles tops and caps.

↳ A large plate mounted with rubberised fingers can accept samples from 10 mm to 240 mm.

↳ It is possible to measure the torque in both clockwise and anti clockwise sense.

↳ The display is an CENTOR FIRST.

↳ For small sized samples, the MICROTORK is made of a little plate with rubberised fingers able to maintain bottles up to 70 mm diameter.



---

**Technical features:**

- ✓ Capacity: 10 Nm
- ✓ Units: Nm – Kgcm - lbin
- ✓ Resolution: 0.01 Nm
- ✓ Accuracy: 0,5% PE
- ✓ Data acquisition rate: 1000Hz

<b><i>Reference</i></b>	ANDITORK FT10
<b><i>Capacity</i></b>	10 Nm
<b><i>Maxi diameter</i></b>	240 mm