

1

2

3

4

STENTOR II

A	CAPACITY:	1/2/5KN	200/500/1000lb
	TRAVEL:	250/350/350mm	9.84/13.8/13.8in
	WEIGHT:	30/40/45KG	66/88/100lb
	HEIGHT (A):	835/935/935mm	32.9/36.8/36.8 in
	MAXIMUM HEIGHT (B)	400/445/445mm	15.7/17.5/17.5in

DIMENSIONS TABLE

WORKABLE WIDTH:	215x320mm	8.47x12.6in
THREAD ON TABLE:	1x M12/4x M8/4x M10	

B


B

C

C

D

D



E

E

4x M10x1.5 ∇ 15[.59]
on \square 150x150[5.9x5.9]

4x M8x1.25 ∇ 15[.59]
on \square 60x60[2.36x2.36]


1x M12x1.75 on table
12.7mm ∇ THRU for loadcell

112.22
[4.418"]

Available surface on the table
 \square 215x320
[8.47"x12.6"]

F

F

Drawn name: FICHE PRODUIT		Assembly name: STENTOR II	
		Date of creation: 14/04/2016	
		Drawn by: CA	Item: STENTOR II

Rev: 2	Date: 08/11/2016	Description:
-----------	---------------------	--------------

1

2

3

4